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## ENGLISH

---

Thank you for purchasing this **LRV7 laboratory food pH electrode**.

For best results with your electrode:

- **read** these operating instructions carefully,
- **comply** with the precautions for use.

The CE marking indicates conformity with European directives, in particular LVD and EMC.

The rubbish bin with a line through it indicates that, in the European Union, the product must undergo selective disposal in compliance with Directive WEEE 2012/19/UE. This equipment must not be treated as household waste.

## DELIVERY CONDITION

---

The food pH electrode is delivered in a cardboard box with:

- a storage bottle,
- an S7-BNC cable approximately 1 metre long or a S7-tight 8-point DIN cable approximately 1 metre long, depending on the order,
- a shim,
- a bottle of storage solution.

For the accessories and spares, consult our web site:

[www.chauvin-arnoux.com](http://www.chauvin-arnoux.com)

## USE

---


This laboratory electrode is used to measure the pH of liquid, semi-solid, and solid foods.

- Unscrew the protective cap.
- Connect the electrode to the pH-meter using the cable supplied with it, S7-BNC or S7- DIN 8. Screw the S7 connector of the cable to the S7 connector of the electrode.
- Remove the storage bottle.
- Calibrate the electrode with buffer solutions, rinsing with de-ionized water between measurements.
- To make a measurement, dip the end of the electrode into the solution to be measured; make sure that the junction is immersed. With a semi-solid or solid food, insert the tip of the electrode into the food as far as the junction. To avoid damaging the electrode, drill a hole in a solid food before inserting the electrode.
- Wait for the measurement to stabilize.
- Rinse the electrode with de-ionized water or a suitable cleaning solution between measurements.
- When the measurements are over, put the electrode back in the storage bottle.
- Unscrew the connecting cable from the pH-meter and screw the protective cap back on.

To avoid damaging the electrode, do not make measurements in any of the following solutions:

- hydrofluoric acid solution,
- acid and alkaline solutions at high concentrations (> 2 mol/L),
- clogging solution (resin, paint, etc.).

## PRESENTATION

## CHARACTERISTICS

Reference: polymer.  
Junction: ceramic and single pore.  
pH measurement range: 0 to 14.  
Temperature range: 0 to 60°C.

Body of the electrode: PVC  
Dimensions: 132 x Ø 16 mm  
Weight: approximately 30 g  
Connector: S7

## CLEANING

---

The presence of crystallized white deposits on the electrode is a normal consequence of the evaporation of the storage solution. To eliminate a deposit, simply rinse the electrode with de-ionized water, then wipe dry with a clean, lint-free cloth.

Do not use aggressive and/or abrasive products (hydrofluoric acid, acetone, etc.). The cleaning operation must not scratch the glass membrane.

Contamination of the junctions is the commonest cause of measurement problems.

- Remove any impurities deposited on the surface of the glass membrane and on the junctions using de-ionized water, or a soft cloth if the deposit is large.
- The commonest cause of clogging is proteins (dairy products, meat). If rinsing is not enough, use a solution of pepsin + 0.1 mol/L HCl.

## MAINTENANCE

---

Keep your electrode in its bottle in a clean 3 mol/L KCl electrolytic solution. Change it as soon as it becomes cloudy.

Do not keep an electrode dry or in distilled or de-ionized water. The membrane may be damaged and the life of the electrode shortened.

Keep the connector clean and dry. If necessary, clean it using a soft cloth wet with soapy water. Rinse using a damp cloth and dry quickly with a dry cloth or forced air.

Regularly check the condition of your electrode and the cleanliness of the junction and of the glass membrane.

## WARRANTY

---

Electrodes are wear parts. The life of your electrode depends on how you use and maintain it.

Except as otherwise stated, our warranty is valid for **12 months** starting from the date on which the equipment was sold. Extract from our General Conditions of Sale, provided on request.

The warranty does not apply in the following cases:

- traces of corrosion on the inside of the connector;
- traces of clogging agents (glue, paint, resin, etc.) on the body, the glass ball or the junctions;
- breakage of the glass;
- damage to the body of the electrode (deformation, discoloration).

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